

# PetriSlides™ for contamination analysis



Substitute for Petri Dishes. Holds filter securely in place. Transparent cover allows microscopic examination without removal. Rectangular base has rounded corners for mounting on microscope stage.

Sterilized storage and examination containers for culturing microorganisms on 47 mm membranes, absorbent pad included.

## Product Specification

### Materials of Construction

Molded polystyrene dish; pure cellulose pad

### Dimension

51 x 76 mm, base

## Applications

- Protect the integrity of sample membrane during microscopic examination
- Store and protect membrane filters for reference

## Product

- ✓ PD1504700 PetriSlides™ for Contamination Analysis 100/pk.

