

COMPACT WHITE LED RING ILLUMINATOR

Model : LCDR48E

- 48 pcs. of ultra bright white LEDs.
- Instant startup
- Flickering Free
- Long life span
 - More than 50,000 hours
- Low power consumption
 - power consumption <4W
- Low heat generation
 - No risk of injury or damage to work pieces by heat radiation
- Dimmable
 - Can be dimmed over a wide range.
- Compact size
 - Single wire direct to wall plug with 12vdc power supply adapter
- Warranty 1 year

Specification:

Electrical

Input voltage : AC 100 ~ 240V

Power consumption : Max <4W

Brightness : ~23,000 lux (height 100mm)



Contact

**Green
Innotec**

Green Innotec Co., Ltd.

99/140 Moo11, Rachadewa, Bangplee

Samutprakarn 10540 Thailand

Tel: +66 (0) 2178 2553


Fax: +66 (0) 25326102

Email: info@greeninnotec.com

www.greeninnotec.com

Mechanical Specification.

Inner diameter : \varnothing 60.5mm
Outer diameter : \varnothing 76mm
Height : 27mm



Contact

**Green
Innotec**

Green Innotec Co., Ltd.

99/140 Moo11, Rachadewa, Bangplee

Samutprakarn 10540 Thailand

Tel: +66 (0) 2178 2553

Fax: +66 (0) 25326102

Email: info@greeninnotec.com

www.greeninnotec.com